CRICOTHYROIDOTOMY TRAINER REVIEW

The DAS 2015 guidelines for the management of unanticipated difficult intubation in adults advocate scalpel cricothyroidotomy for emergency front-of-neck access in the anaesthetised and paralysed patient where intubation and oxygenation has failed (Plan D).

Every anaesthetist should be trained to perform the Plan D technique. This document provides a summary of the features and fidelity of simulator trainers currently available to help provide this training.



We advocate Plan D part task training in two phases: Initial technical and motor skill development with low fidelity simulators, followed by simulation to advance non-technical and psychological skills in a stressful situation. This can be achieved using commercially available simulation apps (e.g., SimMon) in conjuction with the trainers below.

HIGH-FIDELITY EXTREME AIRWAY TRAINER



- Cost £ Set-up Fidelity + CTM Palpable ✓ Vertical Incision ✓ Ultrasound ✓ Latex X
- *Incorporates home-made element* manufacture and assembly. Excellent haptic feedback from the ADAMgel, allows cut-down technique. Only ultrasound-compatible model in this list.

VARADAY ET AL HOMEMADE CRICOTHYROTOMY MODEL

- Set-up Fidelity _ CTM Palpable 🗸 Vertical Incision X Ultrasound X Latex **x**
- Requires assembly, from readily available equipment. Would make ideal model for acquiring initial motor skills for Plan D technique before moving onto higher fidelity models.

SIMULAIDS CRICOTHYROTOMY TRAINER DELUXE No. 135



- Cost **£** £ Set-up Fidelity -CTM Palpable ✓ Vertical Incision X Ultrasound X Latex \(\strice{1} \)
- Internal tube simulates CTM but when present has tendency to redirect bougie into false passage. Very prominent cricothyroid gap. Lubrication required.



- Set-up Fidelity -CTM Palpable 🗸 Vertical Incision X Ultrasound X Latex 🗸
- Very prominent thyroid cartilage with prominent cricothyroid gap. Low-cost foam skin inserts.

BM CRICO-TRAINER FROVA



- Set-up Fidelity -CTM Palpable X Vertical Incision X Ultrasound X Latex ./
- Improved laryngeal cartilage design, low-cost foam skin inserts with very easy set-up.



AERDAL CRICOID STICK SIMULATOR Cost **fff** Unable to perform Plan D tech-Set-up nique on this model with its current Fidelity desian. CTM Palpable X Vertical Incision X Ultrasound X Latex ✓

PHARMABOTICS CRICOTRACHEOTOMY TRAINER



Cost **£££** Set-up Fidelity +

CTM Palpable ✓ Vertical Incision X Ultrasound X

Convincing differentiation of 'cartilaginous' structures and CTM. Lubrication required (not supplied).



eNASCO LIFE/FORM CRICOTHYROTOMY SIMULATOR

Cost **££££** Set-up

Fidelity + CTM Palpable 🗸

Vertical Incision XUltrasound X

Latex X

Only the 'Open' style laryngeal insert permits Plan D technique. Internal clear tube simulates CTM, but when present has tendency to redirect bougie into false passage. Lubrication required (supplied).

SMITHS MEDICAL TRACHEOSTOMY HEAD

CTM Palpable X

Ultrasound x

Latex 🗸

Vertical Incision X



Set-up

Latex OTBC

tion required (not supplied). Fidelity +

Cost **ffff** Rigid thyroid cartilage. Lubrica-



Cost **££££**

Set-up Fidelity +

CTM Palpable 🗸 Vertical Incision X

Ultrasound X Latex 🗸 Lower fidelity than alternative TruCorp product. Laryngeal inserts provide realistic CTM haptic feedback.

VBM CRICO-TRAINER PIG



Cost **ffff** Set-up

Fidelity + CTM Palpable 🗸

Vertical Incision X Ultrasound X

Latex /

Requires separate supply of porcine larynx, with ensuing refrigerated storage and hygiene issues.

AERDAL DELUXE DIFFICULT AIRWAY TRAINER



Fidelity +

CTM Palpable 🗸 Vertical Incision X

Ultrasound X Latex 🗸

Cost **££££** Optional CTM tape available Set-up OTBC to simulate fibrous membrane for puncture. Allows combined airway management and Plan D simulation.

SYNDAVER ADULT CRIC TRAINER



Fidelity + CTM Palpable ✓ Vertical Incision ✓

Ultrasound X

Latex X

Cost **ffff** Not tested directly; requires Set-up OTBC shipping direct from supplier in United States, 'Low-cost Adult Cric Trainer' no longer available from this manufacturer.

TRUCORP AIRSIM ADVANCE COMBO



Set-up Fidelity + Vertical Incision X

Ultrasound X

Latex 🗸

Cost **ffff** Floating larynx provides additional realism. Allows combined advanced airway management and Plan D simulation. Lubrication required (supplied).

KEY

● INVOLVED ● STRAIGHTFORWARD ● VERY EASY

- LOW: skin, thyroid cartilage, cricoid cartilage, hollow trachea

+ HIGH: as above incorporating: skin texture, subcutaneous layers, fibrous CTM, tracheal rings, mandible, face, tracheal path